



AibGenesis™ Mouse Anti-VACV C9L Monoclonal Antibody

Cat. No.: VRS-0224-YT313

This product is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Overview

Target	C9L
Specificity	VACV
Clone	20749YF
Host Species	Mouse
Antibody Isotype	IgG
Species Reactivity	Vaccinia virus

Product Properties

Immunogen	Synthetic peptide sequence
Purity	>95%
Purification	Protein A or G affinity chromatography
Concentration	Lot specific

Packaging, Storage & Formulations

Form	Liquid, please specify if lyophilized powder is required.
Formulation	PBS (pH 7.4) without any preservative, please specify if any additional requirements.
Preservative	Preservative-free, please specify if any additional requirements.
Storage	Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles. Refer to the COA file for specifics.

Applications

Application	ELISA
-------------	-------

Application Notes	The antibody is screened by ELISA using the synthetic peptide corresponding to a partial sequence of Vaccinia virus C9L as the capture antigen. Additional applications have not been validated. It is recommended that the end users validate other applications and determine optimal concentrations/dilution under their experimental conditions.
-------------------	--

Other Product Details

Type	Primary Antibody
Clonality	Monoclonal
ProAb	Made to order

Virus Details

Virus Classification	Double-stranded DNA virus
Virus Family	<i>Poxviridae</i>
Virus Subfamily	<i>Chordopoxvirinae</i>
Virus Genus	<i>Orthopoxvirus</i>
Species	<i>Vaccinia virus</i>
Accession number	AY 243312
RefSeq number	NC_6998
Virus Abbrev	VACV
Virus Name	Vaccinia virus
Virus Alternatives Names	Vaccinia Virus; VACV
Genome Composition	dsDNA

Target

Target Alternative Names	VACWR019, Vaccinia virus, VACV, C9L
Gene ID	3707634
UniProt ID	P17372



AibGenesis™

is an advanced AI-driven platform designed to create novel antibody sequences with unprecedented speed and precision. By integrating deep learning, structure prediction, and comprehensive immunological datasets, **AibGenesis™** intelligently designs antibodies optimized for affinity, stability, and developability. The platform generates antibody products that support basic scientific research, drug development, and diagnostic applications.