

## AibGenesis™ ViroAb™ Mouse Anti-Human CD147 Monoclonal Antibody (XS-0195)

Cat. No.: VRS-0224-YT777

This product is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Overview

Target	CD147
Specificity	This antibody recognizes a 50-60 kDa cell surface protein which is the Human homologue of Rat OX-47 antigen. The antibody recognizes the 118 kDa cell surface glycoprotein CD146, also known as MUC18, Mel-CAM and S-endo.
Clone	XS-0195
Host Species	Mouse
Antibody Isotype	IgG1
Species Reactivity	Human

### Product Properties

Immunogen	Recombinant full length protein (Human)
Purity	>95% as determined by SDS-PAGE
Purification	Antigen affinity chromatography
Concentration	Lot specific

### Packaging, Storage & Formulations

Form	Liquid
Formulation	PBS (pH: 7.2), 1% BSA, 0.09% Sodium Azide
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. Do not freeze. Refer to the COA file for specifics.

### Applications

Application	FC
-------------	----

Application Notes 20 µl reagent for 10<sup>6</sup> cells in a suspension. The optimal working dilutions should be determined by the end user.

## Other Product Details

Clonality	Monoclonal
Related Disease	Malignant Tumors, Infectious Diseases, Rheumatoid Arthritis, Myocardial Infarction, Multiple Sclerosis, Renal Fibrosis
Disease Classification	Immune System Diseases

## Virus Details

Virus Classification	Positive-sense single-stranded RNA Virus
Virus Family	<i>Coronaviridae</i>
Virus Subfamily	<i>Orthocoronavirinae</i>
Virus Genus	<i>Betacoronavirus</i>
Species	<i>Severe acute respiratory syndrome-related coronavirus</i>
Accession number	<a href="#">AY 274119</a>
RefSeq number	<a href="#">NC 4718</a>
Virus Abbrev	SARS-CoV
Virus Name	Severe acute respiratory syndrome coronavirus
Virus Alternatives Names	Severe acute respiratory syndrome-related coronavirus; SARS-CoV; SARS coronavirus; Severe acute respiratory syndrome-related coronavirus 2; SARS-CoV-2
Genome Composition	ssRNA(+)

## Target

Introduction CD147, also known as neurothelin or basigin, is a member of the Ig superfamily. It is a 55-65 kD type I transmembrane glycoprotein which is primarily expressed on leukocytes, erythrocytes, platelets, and endothelial cells. CD147 is reported to have a function during embryonal brain development and/or play a role in integrin-mediated adhesion in brain endothelia.

Target Alternative Names Basigin, BSG, SLC7A11, 5F7, OK blood group antigen, OK, Collagenase stimulatory factor, Extracellular matrix metalloproteinase inducer, EMMPRIN, Leukocyte activati

on antigen M6, M6, Tumor cell-derived collagenase stimulatory factor, TCSF, CD147 antigen, CD147


---

Gene ID [682](#)

---

UniProt ID [P35613](#)

---



**AibGenesis™**

is an advanced AI-driven platform designed to create novel antibody sequences with unprecedented speed and precision. By integrating deep learning, structure prediction, and comprehensive immunological datasets, **AibGenesis™** intelligently designs antibodies optimized for affinity, stability, and developability. The platform generates antibody products that support basic scientific research, drug development, and diagnostic applications.